4-Person Team Division - Medical Imaging Hackathon

Team Structure & Responsibilities

Person 1: AWS Infrastructure & Backend Lead

Role: Cloud Infrastructure & API Development Time Allocation: Full 3 days

Person 2: AI/ML & SageMaker Specialist

Role: Machine Learning Models & Al Integration Time Allocation: Full 3 days

Person 3: IBM watsonx.ai & LLM Integration

Role: Large Language Model & Clinical Report Generation Time Allocation: Full 3 days

Person 4: Frontend & Integration Lead

Role: User Interface & System Integration Time Allocation: Full 3 days

PERSON 1: AWS Infrastructure & Backend Lead

Day 1 Tasks (Setup & Storage)

Morning (9 AM - 12 PM)

AWS Account & Basic Setup

bash

1. Create AWS Account and configure CLI

aws configure

aws sts get-caller-identity # Verify setup

2. Create S3 bucket for medical images

aws s3 mb s3://hackathon-medical-images-\$(date +%s)

aws s3api put-bucket-versioning --bucket your-bucket --versioning-configuration Status=Enabled

3. Enable S3 event notifications

aws s3api put-bucket-notification-configuration \


- --bucket your-bucket \
- --notification-configuration file://s3-notification.json

Create s3-notification.json:

json

Lambda Functions Development

Main Processing Lambda (medical-processor):



```
import json
import boto3
import base64
import os
from datetime import datetime
s3_client = boto3.client('s3')
lambda_client = boto3.client('lambda')
def lambda_handler(event, context):
  try:
    # Extract S3 event information
    bucket = event['Records'][0]['s3']['bucket']['name']
    key = event['Records'][0]['s3']['object']['key']
    print(f"Processing image: {key} from bucket: {bucket}")
    # Get image metadata
    response = s3_client.head_object(Bucket=bucket, Key=key)
    image_size = response['ContentLength']
    # Create processing job
    job_data = {
      'bucket': bucket,
      'key': key,
      'timestamp': datetime.utcnow().isoformat(),
      'size': image_size,
      'status': 'processing'
    # Store job info in S3
    job_key = f"jobs/{key.split('/')[-1]}.json"
    s3_client.put_object(
       Bucket=bucket,
      Key=job_key,
       Body=json.dumps(job_data),
       ContentType='application/json'
    # Trigger AI analysis (invoke Person 2's function)
    lambda_client.invoke(
       FunctionName='ai-analysis-function',
      InvocationType='Event',
       Payload=json.dumps({
         'bucket': bucket,
         'image_key': key,
```

```
'job_key': job_key
})

return {
    'statusCode': 200,
    'body': json.dumps({
        'message': 'Processing started',
        'job_id': job_key
    })
}

except Exception as e:
    print(f"Error: {str(e)}")
    return {
        'statusCode': 500,
        'body': json.dumps({'error': str(e)})
}
```

Deploy Lambda:

```
bash

# Create deployment package
zip -r medical-processor.zip lambda_function.py

# Create Lambda function
aws lambda create-function \
--function-name medical-processor \
--runtime python3.9 \
--role arn:aws:iam::account:role/lambda-execution-role \
--handler lambda_function.lambda_handler \
--zip-file fileb://medical-processor.zip
```

Day 2 Tasks (API Gateway & Integration)

Morning (9 AM - 12 PM)

API Gateway Setup

bash

# Create REST API aws apigateway create-rest-apiname medical-imaging-api	
# Create resources and methods # /upload endpoint for image uploads	
# /status/{job_id} for checking processing status	
# /results/{job_id} for getting final results	

API Gateway Lambda Integration:

(`
	python	

```
# upload-api-handler.py
import json
import boto3
import uuid
import base64
s3_client = boto3.client('s3')
def lambda_handler(event, context):
  try:
    # Handle CORS
    headers = {
      'Access-Control-Allow-Origin': '*',
      'Access-Control-Allow-Headers': 'Content-Type',
      'Access-Control-Allow-Methods': 'POST, GET, OPTIONS'
    if event['httpMethod'] == 'OPTIONS':
      return {
         'statusCode': 200,
         'headers': headers,
         'body': ''
      }
    # Handle file upload
    if event['httpMethod'] == 'POST':
      # Parse multipart form data or base64
      body = json.loads(event['body'])
      # Generate unique filename
      file_id = str(uuid.uuid4())
      filename = f"uploads/{file_id}.jpg"
      # Upload to S3
      s3_client.put_object(
         Bucket=os.environ['BUCKET_NAME'],
         Key=filename,
         Body=base64.b64decode(body['image']),
         ContentType='image/jpeg'
      )
      return {
         'statusCode': 200,
         'headers': headers,
         'body': json.dumps({
           'job_id': file_id,
```

```
'status': 'uploaded',
         'message': 'Image uploaded successfully'
      })
  # Handle status check
  elif event['httpMethod'] == 'GET' and 'job_id' in event['pathParameters']:
    job_id = event['pathParameters']['job_id']
    # Check job status in S3
      response = s3_client.get_object(
         Bucket=os.environ['BUCKET_NAME'],
         Key=f"results/{job_id}.json"
      result = json.loads(response['Body'].read())
      return {
         'statusCode': 200,
         'headers': headers,
         'body': json.dumps(result)
      }
    except:
      return {
         'statusCode': 202,
         'headers': headers,
         'body': json.dumps({'status': 'processing'})
      }
except Exception as e:
  return {
    'statusCode': 500,
    'headers': headers,
    'body': json.dumps({'error': str(e)})
```

Monitoring & Error Handling

CloudWatch Setup:

python

Day 3 Tasks (Final Integration & Testing)

Full Day (9 AM - 6 PM)

System Integration & Testing

- 1. Connect all Lambda functions
- 2. Test end-to-end flow
- 3. Performance optimization
- 4. Error handling refinement
- 5. Documentation for team

Final API Endpoints Documentation:

```
POST /upload - Upload medical image
GET /status/{job_id} - Check processing status
GET /results/{job_id} - Get analysis results
```

PERSON 2: AI/ML & SageMaker Specialist

Day 1 Tasks (Model Research & Setup)

Morning (9 AM - 12 PM)

SageMaker Environment Setup

```
python

# setup-sagemaker.py
import boto3
import sagemaker
from sagemaker import get_execution_role

# Initialize SageMaker session
sagemaker_session = sagemaker.Session()
role = get_execution_role() # Create this IAM role

# Create S3 bucket for models
bucket = sagemaker_session.default_bucket()
prefix = 'medical-imaging-models'

print(f"SageMaker role: {role}")
print(f"S3 bucket: {bucket}")
```

Create SageMaker Execution Role:

```
# IAM role for SageMaker

aws iam create-role --role-name SageMakerExecutionRole \
--assume-role-policy-document file://trust-policy.json

aws iam attach-role-policy --role-name SageMakerExecutionRole \
--policy-arn arn:aws:iam::aws:policy/AmazonSageMakerFullAccess
```

Afternoon (1 PM - 6 PM)


Model Selection & Preparation

Option 1: Use Pre-trained HuggingFace Model

python			

```
# medical-model-setup.py
from sagemaker.huggingface import HuggingFaceModel
import sagemaker
def deploy_medical_model():
  # Use a pre-trained medical imaging model
  huggingface_model = HuggingFaceModel(
    model_data="s3://huggingface-models/medical-imaging/",
    role=get_execution_role(),
    transformers_version="4.21",
    pytorch_version="1.12",
    py_version="py39",
    predictor_cls=sagemaker.predictor.Predictor
  # Deploy to endpoint
  predictor = huggingface_model.deploy(
    initial_instance_count=1,
    instance_type="ml.m5.large",
    endpoint_name="medical-imaging-endpoint"
  return predictor
```

Option 2: Quick Custom Model (Faster for hackathon)



```
# simple-medical-classifier.py
import json
import torch
import torchvision.transforms as transforms
from PIL import Image
import boto3
class SimpleMedicalClassifier:
  def ___init___(self):
    # Load pre-trained ResNet and adapt for medical imaging
    self.model = torch.hub.load('pytorch/vision:v0.10.0', 'resnet50', pretrained=True)
    self.model.eval()
    self.transform = transforms.Compose([
       transforms.Resize(256),
      transforms.CenterCrop(224),
      transforms.ToTensor(),
      transforms.Normalize(mean=[0.485, 0.456, 0.406],
                  std=[0.229, 0.224, 0.225])
    ])
    # Medical conditions mapping (simplified for demo)
    self.conditions = {
       0: "normal",
      1: "potential_abnormality",
       2: "urgent_finding"
  def analyze_image(self, image_bytes):
    # Convert bytes to PIL Image
    image = Image.open(io.BytesIO(image_bytes))
    # Apply transforms
    input_tensor = self.transform(image).unsqueeze(0)
    # Get prediction
    with torch.no_grad():
       outputs = self.model(input_tensor)
       probabilities = torch.nn.functional.softmax(outputs[0], dim=0)
    # Return structured results
    return {
       "findings": self.conditions.get(torch.argmax(probabilities).item(), "unknown"),
       "confidence": float(torch.max(probabilities)),
       "all_probabilities": probabilities.tolist()[:3] # Top 3
```

```
# Lambda function for AI analysis
def lambda_handler(event, context):
  s3_client = boto3.client('s3')
  classifier = SimpleMedicalClassifier()
  try:
    bucket = event['bucket']
    image_key = event['image_key']
    job_key = event['job_key']
    # Download image from S3
    response = s3_client.get_object(Bucket=bucket, Key=image_key)
    image_bytes = response['Body'].read()
    # Analyze image
    results = classifier.analyze_image(image_bytes)
    # Add metadata
    results['image_key'] = image_key
    results['processing_time'] = '2.3 seconds' # Mock for demo
    results['model_version'] = 'v1.0'
    # Save results to S3
    result_key = f"ai-results/{job_key}"
    s3_client.put_object(
      Bucket=bucket,
      Key=result_key,
      Body=json.dumps(results),
      ContentType='application/json'
    # Trigger LLM processing (Person 3's function)
    lambda_client = boto3.client('lambda')
    lambda_client.invoke(
      FunctionName='Ilm-processing-function',
      InvocationType='Event',
      Payload=json.dumps({
        'bucket': bucket,
        'ai_results_key': result_key,
        'job_key': job_key
      })
    return {
      'statusCode': 200,
      'body': json.dumps(results)
```

```
except Exception as e:
    print(f"Al Analysis Error: {str(e)}")
    return {
        'statusCode': 500,
        'body': json.dumps({'error': str(e)})
    }
```

Day 2 Tasks (Model Deployment & Testing)

Morning (9 AM - 12 PM)

Deploy & Test Models

```
bash


# Deploy Lambda function for AI processing

zip -r ai-analysis.zip lambda_function.py model_files/

aws lambda create-function \
--function-name ai-analysis-function \
--runtime python3.9 \
--role arn:aws:iam::account:role/lambda-execution-role \
--handler lambda_function.lambda_handler \
--zip-file fileb://ai-analysis.zip \
--timeout 300 \
--memory-size 1024
```

Afternoon (1 PM - 6 PM)

Model Optimization & Validation



```
# model-validation.py
import json
import time
def validate_model_performance():
  test_cases = [
    "sample_ct_normal.jpg",
    "sample_ct_stroke.jpg",
    "sample_mri_hemorrhage.jpg"
  1
  results = []
  for test_image in test_cases:
    start_time = time.time()
    # Test your model
    result = analyze_test_image(test_image)
    processing_time = time.time() - start_time
    results.append({
      'image': test_image,
      'result': result,
      'processing_time': processing_time
    })
  return results
def create_demo_dataset():
  """Create sample medical images with known results for demo"""
  sample_results = {
    "normal_ct": {
       "findings": "normal",
       "confidence": 0.92,
       "description": "No acute abnormalities detected"
    },
    "stroke_ct": {
      "findings": "urgent_finding",
      "confidence": 0.87,
      "description": "Possible acute stroke - left MCA territory"
    },
    "hemorrhage_ct": {
       "findings": "urgent_finding",
       "confidence": 0.94,
       "description": "Intracranial hemorrhage detected"
```

```
}
return sample_results
```

Day 3 Tasks (Integration & Performance)

Full Day (9 AM - 6 PM)

Final Model Integration & Testing

- 1. Integration testing with Person 1's infrastructure
- 2. Performance optimization
- 3. Mock realistic medical results for demo
- 4. Coordinate with Person 3 for LLM handoff

PERSON 3: IBM watsonx.ai & LLM Integration

Day 1 Tasks (IBM Setup & LLM Access)

Morning (9 AM - 12 PM)

IBM Cloud Setup


```
# Install IBM Cloud CLI
curl -fsSL https://clis.cloud.ibm.com/install/linux | sh

# Login and setup
ibmcloud login --apikey YOUR_API_KEY
ibmcloud target -r us-south -g default
```

Python Setup for watsonx.ai:

python			
python			

```
# watson-setup.py
from ibm_watson_machine_learning import APIClient
import json
# Watson ML credentials
wml_credentials = {
  "url": "https://us-south.ml.cloud.ibm.com",
  "apikey": "your-api-key-here",
  "instance_id": "your-instance-id"
# Initialize client
client = APIClient(wml_credentials)
# Set project
client.set.default_project("your-project-id")
# List available foundation models
models = client.foundation_models.get_model_specs()
print("Available models:", models)
# Test connection
print("Watson ML client initialized successfully!")
```



```
# medical-prompt-templates.py
class MedicalPromptTemplates:
  @staticmethod
  def clinical_summary_prompt(ai_findings, patient_context="emergency"):
    return f"""
You are an expert radiologist Al assistant. Based on the medical imaging analysis results, generate a concise
IMAGING FINDINGS:
{ai_findings}
CLINICAL CONTEXT: {patient_context}
Please provide:
1. KEY FINDINGS: (2-3 sentences max)
2. CLINICAL SIGNIFICANCE: (Critical/Moderate/Low concern)
3. RECOMMENDED ACTIONS: (Immediate steps)
4. FOLLOW-UP: (If needed)
Keep the language clear and actionable for emergency physicians. Focus on time-sensitive information.
CLINICAL SUMMARY:
0.00
  @staticmethod
  def treatment_protocol_prompt(findings, condition):
    return f"""
Based on the medical imaging findings showing {condition}, provide the appropriate emergency treatment pr
FINDINGS: {findings}
Provide a structured treatment protocol including:
- IMMEDIATE ACTIONS (within 15 minutes)
- DIAGNOSTIC WORKUP (additional tests needed)
- TREATMENT OPTIONS (evidence-based)
- CONSULTATION REQUIREMENTS (specialist referrals)
TREATMENT PROTOCOL:
0.00
  @staticmethod
  def patient_explanation_prompt(technical_findings):
    return f"""
Convert the following technical medical findings into simple language for patient/family explanation:
TECHNICAL FINDINGS: {technical_findings}
```

```
Provide a clear, compassionate explanation that:
- Uses simple, non-technical language
- Explains what was found
- Explains next steps
- Addresses likely concerns
PATIENT EXPLANATION:
0.00
# Test the prompts
def test_prompts():
  sample_findings = {
    "findings": "urgent_finding",
    "confidence": 0.87,
    "description": "Possible acute stroke - left MCA territory"
  prompt = MedicalPromptTemplates.clinical_summary_prompt(
    json.dumps(sample_findings),
    "emergency"
  print("Generated Prompt:")
  print(prompt)
```

Day 2 Tasks (LLM Integration & Processing)

Morning (9 AM - 12 PM)

Granite LLM Integration

python

```
# Ilm-processor.py
from ibm_watson_machine_learning import APIClient
import ison
import boto3
class GraniteMedicalLLM:
  def __init__(self):
    self.wml_credentials = {
      "url": "https://us-south.ml.cloud.ibm.com",
      "apikey": os.environ['IBM_API_KEY'],
      "instance_id": os.environ['IBM_INSTANCE_ID']
    self.client = APIClient(self.wml_credentials)
    self.client.set.default_project(os.environ['IBM_PROJECT_ID'])
    # Model parameters
    self.generation_params = {
      "max_new_tokens": 300.
      "temperature": 0.3,
      "top_p": 0.9,
      "repetition_penalty": 1.1
  def generate_clinical_summary(self, ai_findings):
    """Generate clinical summary using Granite LLM"""
    try:
       prompt = MedicalPromptTemplates.clinical_summary_prompt(ai_findings)
      # Generate response
      response = self.client.foundation_models.generate_text(
         model_id="ibm/granite-13b-chat-v2", # or latest Granite model
         prompt=prompt,
         params=self.generation_params
      return {
         "clinical_summary": response,
         "model_used": "granite-13b-chat-v2",
         "confidence": "high",
         "generated_at": datetime.utcnow().isoformat()
      }
    except Exception as e:
      print(f"LLM Generation Error: {str(e)}")
      # Fallback to template-based response
```

```
return self.fallback_summary(ai_findings)
  def fallback_summary(self, ai_findings):
    """Fallback summary in case LLM fails"""
    findings_data = json.loads(ai_findings) if isinstance(ai_findings, str) else ai_findings
    if findings_data.get('findings') == 'urgent_finding':
      return {
         "clinical_summarv": """
KEY FINDINGS: Potential urgent abnormality detected on imaging with high confidence.
CLINICAL SIGNIFICANCE: HIGH CONCERN - requires immediate physician review.
RECOMMENDED ACTIONS:
- Immediate physician evaluation
- Consider specialist consultation
- Monitor patient closely
FOLLOW-UP: Formal radiologist interpretation recommended within 1 hour.
         "model_used": "fallback_template",
         "confidence": "template_based"
      }
    else:
      return {
         "clinical_summarv": """
KEY FINDINGS: No acute abnormalities detected on initial AI screening.
CLINICAL SIGNIFICANCE: LOW CONCERN - routine findings.
RECOMMENDED ACTIONS:
- Continue standard clinical evaluation
- Formal radiologist review as scheduled
FOLLOW-UP: Standard radiology reporting timeframe.
         0.010
         "model_used": "fallback_template",
         "confidence": "template_based"
# Lambda function for LLM processing
def lambda_handler(event, context):
  s3_client = boto3.client('s3')
  Ilm_processor = GraniteMedicalLLM()
  try:
    bucket = event['bucket']
    ai_results_key = event['ai_results_key']
    job_key = event['job_key']
    # Get AI analysis results from S3
    response = s3_client.get_object(Bucket=bucket, Key=ai_results_key)
    ai_results = json.loads(response['Body'].read())
```

```
# Generate clinical summary
  clinical_summary = Ilm_processor.generate_clinical_summary(ai_results)
  # Combine results
  final_results = {
    "ai_analysis": ai_results,
    "clinical_summary": clinical_summary,
    "processing_complete": True,
    "total_processing_time": "2.8 seconds", # Mock for demo
    "timestamp": datetime.utcnow().isoformat()
  # Save final results to S3
  final_key = f"results/{job_key.replace('jobs/', '').replace('.json', '')}.json"
  s3_client.put_object(
    Bucket=bucket,
    Key=final_key,
    Body=json.dumps(final_results),
    ContentType='application/json'
  return {
    'statusCode': 200,
    'body': json.dumps({
      'message': 'Clinical summary generated',
      'results_key': final_key
    })
except Exception as e:
  print(f"LLM Processing Error: {str(e)}")
  return {
    'statusCode': 500,
    'body': json.dumps({'error': str(e)})
```

Medical Knowledge Base

python

```
# medical-guidelines.py
class MedicalGuidelinesDB:
  def __init__(self):
    self.guidelines = {
       "stroke": {
         "protocol": "Activate stroke code immediately",
         "time_window": "4.5 hours for thrombolysis".
         "imaging": "Non-contrast CT first, then CT angiography",
         "treatment": "Consider tPA if within window"
      },
      "hemorrhage": {
         "protocol": "Neurosurgery consultation stat",
         "imaging": "Non-contrast CT diagnostic",
         "treatment": "Reverse anticoagulation if present",
         "monitoring": "Neuro checks q15min"
      },
      "normal": {
         "protocol": "Standard emergency evaluation",
         "next_steps": "Clinical correlation recommended".
         "follow_up": "Routine radiology review"
      }
  def get_guideline(self, condition):
    return self.guidelines.get(condition.lower(), self.guidelines["normal"])
  def enhance_summary_with_guidelines(self, summary, condition):
    guideline = self.get_guideline(condition)
    enhanced = f"""
{summary}
CLINICAL GUIDELINES:
Protocol: {guideline.get('protocol', 'Standard care')}
Timing: {guideline.get('time_window', 'No specific time constraints')}
Next Steps: {guideline.get('treatment', 'Continue clinical evaluation')}
    return enhanced
```

Day 3 Tasks (Testing & Optimization)

Full Day (9 AM - 6 PM)

LLM Testing & Integration

Create fallback mechanisms Integration testing with team					
PERSON 4: Frontend & Integration Lead					
Day 1 Tasks (Frontend Setup)					
Morning (9 AM - 12 PM)					
Basic HTML/CSS/JS Setup					
html					

1. Test different prompt variations

2. Optimize response times

```
<!-- index.html -->
<!DOCTYPE html>
<html lang="en">
<head>
  <meta charset="UTF-8">
  <meta name="viewport" content="width=device-width, initial-scale=1.0">
  <title>Emergency Medical Imaging Al</title>
  <link href="https://cdnjs.cloudflare.com/ajax/libs/tailwindcss/2.2.19/tailwind.min.css" rel="stylesheet">
  <style>
    .upload-area {
      border: 2px dashed #cbd5e0;
      transition: all 0.3s ease:
    .upload-area.dragover {
      border-color: #3182ce;
      background-color: #ebf8ff;
    .processing {
      animation: pulse 2s infinite:
  </style>
</head>
<body class="bg-gray-100">
  <div class="container mx-auto px-4 py-8">
    <h1 class="text-4xl font-bold text-center text-gray-800 mb-8">
      Emergency Medical Imaging Al
    </h1>
    <!-- Upload Section -->
    <div class="max-w-2xl mx-auto">
      <div id="uploadArea" class="upload-area rounded-lg p-8 text-center mb-6">
        <div id="uploadContent">
          <svg class="mx-auto h-12 w-12 text-gray-400 mb-4" stroke="currentColor" fill="none" viewBox</pre>
            <path d="M28 8H12a4 4 0 00-4 4v20m32-12v8m0 0v8a4 4 0 01-4 4H12a4 4 0 01-4-4v-4m3</p>
          </svg>
          Upload CT or MRI Scan
          Drag and drop or click to select
          <input type="file" id="imageInput" accept="image/*,.dcm" class="hidden">
          <button onclick="document.getElementById('imageInput').click()"</pre>
              class="bg-blue-500 hover:bg-blue-600 text-white px-6 py-3 rounded-lg">
            Select Image
          </button>
        </div>
        <!-- Processing State -->
        <div id="processingState" class="hidden">
```

```
<div class="processing">
      <div class="inline-block w-8 h-8 border-4 border-blue-500 border-l-transparent rounded-full</p>
    Analyzing Medical Image...
    Initializing AI analysis
    <div class="bg-gray-200 rounded-full h-2 mt-4">
     <div id="progressBar" class="bg-blue-500 h-2 rounded-full transition-all duration-500" style=</pre>
   </div>
 </div>
</div>
<!-- Results Section -->
<div id="resultsSection" class="hidden">
  <div class="bg-white rounded-lg shadow-lg p-6">
    <h2 class="text-2xl font-bold text-gray-800 mb-4">Analysis Results</h2>
   <!-- AI Findings -->
    <div class="mb-6">
     <h3 class="text-lg font-semibold text-gray-700 mb-2">AI Detection Results</h3>
     <div id="aiFindings" class="bg-gray-50 rounded p-4">
       <!-- Dynamic content -->
```